

RelyX™ ARC Adhesive Resin Cement

For strength and esthetics, it's an easy choice.



Photo courtesy of Robert A. Lowe, D.D.S.

- Recommended for indirect restorations
- Easy to use
- Dispenses uniform, procedure-specific dosages
- Provides easy excess removal

Superior strength. Easy to use.

RelyX™ ARC Adhesive Resin Cement from 3M ESPE is as easy to work with as it is strong. Featuring 3M ESPE's innovative Clicker™ Dispenser, RelyX ARC adhesive resin cement delivers precise dosages with less waste. It bonds indirect restorations securely and adheres to metals, porcelain, ceramic, pre-cured composite, enamel and dentin—yet the excess comes off easily. You can also use it for endodontic post cementation, amalgam bonding and maryland bridge cementation. And better yet, it is compatible with Adper™ Single Bond Adhesive from 3M ESPE.



	RelyX™ ARC Cement	PANAVIA® 21	C&B-Metabond®
Bond Strength (24 Hour-MPa)			
-Porcelain to Dentin	23.7 (6.7)	7.6 (5.4)	20.9 (5.0)
-Porcelain to Enamel	24.4 (5.4)	25.4 (7.1)	20.2 (7.0)
Compressive Strength (24 Hour-MPa)	327.8 (9.9)	292.2 (13.6)	67.9 (1.5)
Diametral Strength (24 Hour-MPa)	59.8 (2.3)	62.1 (1.9)	22.5 (0.3)

RelyX™ ARC Adhesive Resin Cement

What the Experts Say

Dentists around the world are discovering the benefits of RelyX™ ARC Adhesive Resin Cement from 3M ESPE. See what industry experts say about this product, then try it for yourself to experience the unique advantages for you and your patients.

"Consultants rated RelyX ARC adhesive resin cement very good to excellent for its ease of use, consistency and dispensing."

The Dental Advisor, Vol. 16, 1999

"Dispensing system very convenient and fast; good flow characteristics and color; easy clean-up. 87% of 24 evaluators rated product excellent or good and worthy of trial by colleagues."

An independent research institute, July 1999

"Low viscosity, very easy to mix, but is not runny. Clean-up is also very easy for a resin cement."

Reality Now, No. 111, June 1999

"I think that 3M has a real winner here—a terrific resin cement and a great new dispensing system."

Dental Economics, March 1999

Ordering Information



Introductory Kit: 3415

RelyX™ ARC Adhesive Resin Cement Introductory Kit

3415

Kit includes:

- 1 - 4.5 g Paste Shade Universal (3415A3) in Clicker™ Dispenser*
- 1 - 4.5 g Paste Shade Transparent (3415A1) in Clicker™ Dispenser*
- 1 - 6 ml Adper™ Single Bond Adhesive (3411)*
- 1 - 5 ml RelyX™ Ceramic Primer (2721)*
- 2 - 3 ml Scotchbond™ Phosphoric Etchant Syringes (7523)*
- Etchant Tips (7523T)*, Brush Tips/Handles (1919)*
- Fiber Tips (1919F)*, Disposable Mix Wells (3401)*
- Bendable Brushes (3402)*, Mix Pads (1994P)*
- Mixing Sticks (1994S)*, Mini-Sponges (7522S)*
- Includes illustrated technique guides and instructions for use

*Product can be ordered individually.

A Complete Family of Solutions

Easy to choose, easy to use.

3M ESPE offers a complete line of RelyX™ Cements for all your cementation needs.

- RelyX™ Veneer Cement
- RelyX™ Temporary Cements
- RelyX™ Luting Cement
- RelyX™ ARC Adhesive Resin Cement
- RelyX™ Unicem Self-Adhesive Universal Resin Cement

For more information visit our Web site, or call our technical hotline:

Technical Hotline: 1-800-634-2249 www.3MESPE.com

3M ESPE

Dental Products

3M Center
Building 275-SE-03
St. Paul, MN 55144-1000
USA

3M Canada

Post Office Box 5757
London, Ontario N6A 4T1
Canada
1-800-265-1840 ext. 6229



10% post-consumer
waste paper

Printed in USA
© 3M ESPE 2003

3M, ESPE, RelyX, Clicker, Adper and Scotchbond are trademarks of 3M or 3M ESPE AG. PANAVIA is a trademark of Kuraray Co., Ltd. C&B-Metabond is a trademark of Parkell, Inc. Variolink is a trademark of Ivoclar Vivadent, Inc.

70-2009-3458-9